

# Quantrex®

*Ultrasonic Cleaning Systems*



## Powerful systems for effective cleaning.

The Quantrex® Ultrasonic Cleaning System Line provides super-strength cleaning every time. Versatile enough for a variety of cleaning applications. When used with L&R's speciality formulated solutions, the self-contained, compact unit offers efficient trouble-free cleaning. Each Quantrex Machine comes standard with increased power ... Strength you can see as soon as you turn the unit on.



**Q140**

The Quantrex Product Line is constructed of a vinyl clad metal wrapper and stainless steel tank, using only the best materials and workmanship. Each tank and body is carefully inspected for strength and durability. The Quantrex Line is backed by the most comprehensive warranty in the industry - all defects in materials or workmanship are covered for 2-½ years.

Sleek design fits into any setting, from a Dental Laboratory to a Jeweler's Retail Counter.

### **Quantrex® Features:**

- *Increased power, yet cleans quietly!*
- *Perfect for countertop use*
- *Optional Heater*
- *60 Minute Timer*
- *Full line of accessories available for all your cleaning needs*
- *ETL-CAN/US, CE approved*
- *ISO 9001: 2000*
- *Constructed with vinyl-clad steel and stainless steel*
- *Fourteen quality inspection steps for strength and durability*
- *Equipped with a stainless steel drain with a multi-positional outlet for easy removal of solution (excluding the Q90)*
- *Available in many different sizes*

*L&R Ultrasonics ... For Cleaning Excellence that you can depend on!*



# Quantrex®

## Ultrasonic Cleaning System



Q140

All L&R Quantrex Systems are backed by an unprecedented warranty.  
**THE DOUBLE WARRANTY.**

It covers the entire unit for defects in material and workmanship for 2½ years. Additionally, the housing is warranted for a period of 5 years and the tank-to-transducer bonding is covered for 10 years. 2½, 5, 10 ... that's quality ..... that's L&R.

Model	Description	Input Power	Output	Overall Dimensions	Shipping Weight.	Tank Capacity	Tank Internal Dimensions
	<i>All Units include Stainless Steel Cover</i>	<i>Average Watts 117 VAC 60 Hz</i>	<i>Frequency</i>	<i>Inches &amp; Centimeters L x W x H</i>	<i>Pounds &amp; Kilograms Approx</i>	<i>Quarts &amp; Liters</i>	<i>Inches &amp; Centimeters L X W x H</i>
Q90	Timer standard, Heater optional	55 110 w/heat	43 KHz	7.0 x 6.38 x 8.75 17.8 x 16.2 x 22.2	8 3.63	2 1.9	5.88 x 5.38 x 4.0 14.9 x 13.7 x 10.2
Q140	Timer standard, Heater optional	95 150 w/Heat	43 KHz	10.25 x 6.5 x 8.25 26.0 x 16.5 x 21.0	10 4.3	3.38 3.2	9.38 x 5.38 x 4.0 23.8 x 13.7 x 10.2
Q210	Timer, Drain standard Heater, optional	135 335 w/heat	43 KHz	12.75 x 7.0 x 11.0 32.4 x 17.8 x 27.9	14 6.36	6 5.7	11.75 x 6.0 x 6.0 29.8 x 15.2 x 15.2
Q280	Timer, Drain standard Heater, optional	240 440 w/Heat	43 KHz	12.75 x 10.5 x 12.13 32.4 x 26.7 x 30.8	21 9.3	11 10.4	11.75 x 9.38 x 6.0 29.8 x 23.8 x 15.2
Q310	Timer, Drain standard, Heater optional	320 520 w/Heat	43 KHz	16.5 x 10.0 x 12.13 41.9 x 25.4 x 30.8	26 11.7	13 12.3	15.5 x 9.0 x 6.0 39.4 x 22.9 x 15.2
Q360	Timer, Drain standard, Heater optional	240 440 w/Heat	43 KHz	12.75 x 10.5 x 14.25 32.4 x 26.7 x 36.20	22 9.9	14.5 13.6	11.5 x 9.25 x 8.0 29.2 x 23.5 x 20.3
Q650	Timer, Drain standard, Heater optional	415 815 w/Heat	43 KHz	21.75 x 13.75 x 13.5 55.2 x 34.93 x 34.29	33 15	26 24.6	19.75 x 11.75 x 6.5 50.2 x 29.8 x 16.5



**L&R Manufacturing Company**

577 Elm Street, PO Box 607, Kearny, NJ 07032-0607 USA Telephone: 201.991.5330 Fax: 201.991.5870  
[www.LRultrasonics.com](http://www.LRultrasonics.com) [info@LRultrasonics.com](mailto:info@LRultrasonics.com)

