



PC3

Ultrasonic Cleaning System Now with Great New Look.

When you require maximum cleaning, but have minimal space, the PC3 is the perfect ultrasonic system. This compact ultrasonic cleaning machine features solid state circuitry for powerful, quiet cleaning action, yet it fits easily on a counter or table top. The perfect choice for a second auxiliary unit as well. When used with L&R's specially formulated cleaning solutions, you will enjoy maximum cleaning results.



The PC3 features stainless steel housing and stainless steel tank cover. Optional accessories: stainless steel mesh basket and chairside kit which includes 2-100ml beakers complete with positioning cover provides great versatility. L&R has accessories available for all your cleaning requirements.



PC3 Ultrasonic Cleaning System



PC 1172

Model	PC3
Description	Self-contained, compact, high performance unit
Input Power	117 VAC 50/60 Hz Average Watts 0.2 amp / 22 watts
Voltage	117 Volts
Peak Output Frequency	22 watts 55 kHz
Tank Capacity	18 oz 0.53 ltr
External Dimensions	5.375" L x 4.24" W x 5" D 13.7 cm L x 10.8 cm W x 12.7 cm D
Internal Dimensions	4.75" L x 3.38" W x 2.63" D 12.1 cm L x 8.6 cm W x 6.7 cm D
Shipping Weight	4 lbs 1.81 Kg

