



UltraDose®

GERMICIDAL Hospital Grade Ultrasonic Cleaner

COVID-19 Effective*

*Emerging Viral Pathogen Claim

L&R UltraDose Germicidal Hospital Grade Ultrasonic Cleaner Concentrate has demonstrated effectiveness against viruses similar to COVID-19 Coronavirus on hard, non-porous surfaces. Therefore, UltraDose Germicidal Hospital Grade Ultrasonic Cleaner Concentrate can be used against COVID-19 Coronavirus when used in accordance with the directions for use against Feline Calicivirus (surrogate for Norovirus) on hard, non-porous surfaces.

Clean and disinfect in one easy step with L&R's dual action UltraDose® Hospital Grade Ultrasonic Cleaner Concentrate.

The unique blend of an ultrasonic cleaning solution and an E.P.A. registered hospital grade disinfectant come together in this effective product. UltraDose® Germicidal eliminates dangerous microorganisms while simultaneously ultrasonically cleaning and disinfecting instruments, apparatus and equipment. Help maintain the safety of your office and staff and keep your instruments clean and disinfected with L&R's UltraDose® Germicidal.

- E.P.A. registered Product
- For use in an ultrasonic cleaner
- Disinfects instruments and the ultrasonic bath
- Cleans non-porous instrument and appliance surfaces, apparatus and equipment
- Performs as an ultrasonic cleaning solution, disinfectant, deodorant, germicidal, fungicidal, virucidal and tuberculocidal all-in-one
- Concentrated; one ounce makes one gallon. Non-irritating and fragrance free

Don't take infection control for granted!

Available in 16 oz. containers.
Each container makes
16 gallons of solution.



because clean matters.

L&R Manufacturing Company
577 Elm Street, P.O. Box 607
Kearny, NJ 07032-0607 USA
Tel: 201.991.5330 Fax: 201.991.5870


www.LRultrasonics.com info@LRultrasonics.com



UltraDose[®]

GERMICIDAL

Hospital Grade Ultrasonic Cleaner



- ✓ Fragrance free
- ✓ Super concentrated
- ✓ Dual purpose proprietary formulation
- ✓ Promotes office staff safety

UltraDose[®] Germicidal is proven effective against viruses much tougher to inactivate than Novel Coronavirus, such as Feline Calicivirus (surrogate for Norovirus). This means UltraDose[®] Germicidal can be used against COVID-19 Coronavirus on hard, non-porous surfaces.

Active Ingredient: Phenylphenol


Container Size: 16 oz.

Dilution: 1 oz./gallon

Gallons/container: Makes 16 gallons of ready to use solution

E.P.A. Registration #: 34810-31-72869

Product Code: UD036



- ⊘ COVID-19*
- ⊘ HIV-1
- ⊘ Pseudomonas Acrusiosa
- ⊘ Trichophyton
- ⊘ Salmonella
- ⊘ Staphylococcus
- ⊘ Seratia Marcescens
- ⊘ Herpes Simplex 2
- ⊘ Mycobacterium Bovis (TB)

because clean matters.

L&R Manufacturing Company
577 Elm Street, P.O. Box 607
Kearny, NJ 07032-0607 USA
Tel: 201.991.5330 Fax: 201.991.5870

www.LRultrasonics.com info@LRultrasonics.com

ISO 9001:2015

SOL.102.0520C