



Quantrex®

The Quantrex® line of Ultrasonic Cleaning Systems sets the standard for precision cleaning; providing versatility, durability and exacting results. Designed for countertop use, the Quantrex series is quiet, self-contained and offers a variety of sizes for every cleaning application, including instruments, parts and apparatus. A full line of accessories is available for all your cleaning requirements. When combined with L&R's premium quality ultrasonic cleaning solutions, Quantrex brand Ultrasonic Cleaning Systems are the superior choice...because clean matters.

Ultrasonic Cleaning Systems

Powerful, durable machines for effective cleaning

- Quick and powerful cleaning
- Designed for countertop use
- Optional heater available
- Standard 60 minute timer
- Shipped complete with cover
- Equipped with a multi-directional stainless steel drain (excluding the Q90)
- ETL-CAN/US, CE approved
- ISO 9001-2008

Quantrex 310

**NOW IN
STAINLESS
STEEL**

Quantrex 140



because clean matters.

L&R Manufacturing Company
577 Elm Street, P.O. Box 607
Kearny, NJ 07032-0607 USA
Tel: 201.991.5330 Fax: 201.991.5870

www.LRultrasonics.com info@LRultrasonics.com



Quantrex®



Quantrex Q140 with Optional Basket

Ultrasonic Cleaning Systems for powerful, reliable and repeatable results

The Quantrex® product line is constructed with a stainless steel housing and tank, using only the best and most durable materials and workmanship. Each machine goes through a series of tests to ensure quality, making L&R the choice of professionals. Every Quantrex machine is backed by the most comprehensive warranty in the industry; all defects in materials and workmanship are covered for 2½ years.

Model *	Input Power Average Watts 117 VAC 60 Hz	Overall Dimensions (L x W x H)	Tank Capacity	Tank Internal Dimensions (L x W x H)
Quantrex 90 **	55 110 with Heat	7.0 x 6.4 x 8.8 in 17.8 x 16.2 x 22.2 cm	0.5 gal 1.9 L	5.9 x 5.4 x 4.0 in 14.9 x 13.7 x 10.2 cm
Quantrex 140	95 150 with Heat	10.3 x 6.5 x 8.3 in 26.0 x 16.5 x 21.0 cm	0.85 gal 3.2 L	9.4 x 5.4 x 4.0 in 23.8 x 13.7 x 10.2 cm
Quantrex 210	135 335 with Heat	12.8 x 7.0 x 11.0 in 32.4 x 17.8 x 27.9 cm	1.5 gal 5.7 L	11.8 x 6.0 x 6.0 in 29.8 x 15.2 x 15.2 cm
Quantrex 310 (No Heater)	320	16.5 x 10.0 x 12.1 in 41.9 x 25.4 x 30.8 cm	3.25 gal 12.3 L	15.5 x 9.0 x 6.0 in 39.4 x 22.9 x 15.2 cm
Quantrex 360	240 440 with Heat	12.8 x 10.5 x 14.3 in 32.4 x 26.7 x 36.2 cm	3.6 gal 13.6 L	11.5 x 9.3 x 8.0 in 29.2 x 23.5 x 20.3 cm
Quantrex 650	415 815 with Heat	21.8 x 13.8 x 12.7 in 55.2 x 34.9 x 32.2 cm	6.50 gal 24.6 L	19.8 x 11.8 x 6.5 in 50.2 x 29.8 x 16.5 cm

**L&R's Triple Warranty covers the entire unit for defects in material and workmanship for 2½ years.
The housing is warrantied for a period of 5 years and the tank-to-transducer bonding is covered for 10 years.**

Quantrex Systems – Shown with Heat

All systems operate at 43kHz.
All systems are available with
international voltages.



* All models include timer, drain and cover. Heater available on most systems. Baskets and accessories sold separately.
** Quantrex 90 does not include drain.

because clean matters.

L&R Manufacturing Company
577 Elm Street, P.O. Box 607
Kearny, NJ 07032-0607 USA
Tel: 201.991.5330 Fax: 201.991.5870

www.LRultrasonics.com info@LRultrasonics.com



Product code: MAC.100.0517.R